SEFFIHAIR®

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DEVICE DESCRIPTION AND USE

The SEFFIHAIR® kit is designed for **physicians** who want to perform the regenerative therapy for regeneration, reconstruction and repair of integumentary and musculoskeletal tissues.

SEFFIHAIR® is a disposable sterile kit to be used on a single patient.

SEFFIHAIR® is designed to facilitate and standardize the harvesting, preparation and grafting of subcutaneous adipose tissue. Several studies have shown that the subcutaneous adipose tissue is made of adipose cells and stromal tissue (SVF) containing mesenchymal cells and especially adipose derived stem cells (ADSCs).

Under this procedure the tissue harvested and prepared with SEFFIHAIR® is implanted in the same patient (autologous graft) in order to stimulate hair growth, improve the microcirculation (angiogenetic action) and an anti-inflammatory and antifibrotic action.

This device can be used on different body areas on the same patient and within the same procedure.

ABSOLUTE CONTRAINDICATIONS

- Ongoing infections in the area of tissue harvesting or grafting.
- Presence of malignancies in the area of tissue harvesting or grafting.
- Pregnancy or breast-feeding.
- Anticoagulant therapy or severe coagulation disorder.
- Allergy to local anesthetic.
- Dysmorphophobia.
- Ongoing immunosuppressive therapies.
- Debilitated subjects.

RELATIVE CONTRAINDICATIONS

- Strong tobacco addiction.
- Presence of non-resorbable material in harvesting or grafting area.

WARNINGS

- 1. This is a disposable device.** Reusing the SEFFIHAIR® system may cause infections and the transmission of communicable diseases.
- 2. This device is sterile.** If the sterile packaging appears damaged, do not use the product.
- 3. This device cannot be resterilized.** Do not autoclave.
- 4. Do not use the device after the expiry date on the label.**
- Before use and during transport of the device, the **surrounding temperature must range between 5°C and 40°C.** Excessive temperatures may compromise the reliability of the device.
- Utmost caution is necessary when using this device on patients in precarious state of health.
- This device is not meant for removal of adipose tissue in a body remodeling treatment or weight reduction.
- During the procedure do not force the guide rod. In case of breakage of the rod or of other parts of the guide **interrupt immediately** the procedure.
- The harvesting procedure must be performed in body areas where the pinch test is = or > 4 cm.
- In order to guarantee the traceability of the product, please put the sticker with QR code and barcode located on the bottom of SEFFIHAIR® box on the form for the informed consent to be signed by the patient both on the form for the informed consent to be signed by the patient and on the medical record.

CAUTIONS

- The use of this device is limited only to those physicians who have been specially trained.
- The procedure outcomes vary depending on the patient's age, state of health, life habits and surgical site as well as on the physician's experience.
- Results achieved may not be definitive.
- The quantity of tissue harvested must be limited to what is actually needed.
- The surgical field needs to be completely sterile.** Please follow medical standards designed to guarantee the maintenance of sterility during the procedure (e.g. do not touch the non-sterile parts of the patient or operator such as garments, jewellery, glasses, phones, etc.).

UNDESIRE EFFECTS

Undesired effects deriving from the procedure may be: asymmetries, infection, inflammation, hematomas, ecchymoses, skin irregularities, scar healing defects, skin dyschromia, fat necrosis, oil cysts, allergic reactions.

KIT CONTENTS SEFFIHAIR®

THE STERILE BAG INSIDE THE BOX CONTAINS:

- no. 1 guide
- no. 1 1.0mm-hole infiltration and harvesting cannula
- no. 1 multiple syringe decanting stand
- no. 1 syringe plunger lock
- no. 3 10 ml luer lock syringes
- no. 3 3 ml luer lock syringes
- no. 3 test tubes
- no. 3 plugs for luer lock syringes
- no. 2 connectors for luer lock syringes
- no. 1 cylindrical container
- no. 8 needles with different caliber and length
- no. 1 drape
- no. 1 drape with opening
- no. 4 gauze swabs
- no. 1 tray

STEPS IN THE PROCEDURE

The procedure is to be performed in sterile conditions: the skin should be carefully cleaned with a disinfecting solution both in the harvesting area and in the implantation area prior to the procedure. Do not compromise the sterility conditions during the procedure. You should absolutely avoid any contact with non-sterile parts of the patient or operator (garments, jewellery, glasses, phones, etc.)

HARVESTING

The area where the tissue is harvested should be chosen based on the following criteria:

- adequate quantity of subcutaneous adipose tissue;
- size of the area: the harvesting area should be at least 20 cm in diameter;
- the most common harvesting areas are: abdomen, hips and the trochanteric region.

Harvesting phases:

A. PREPARATION

- A1.** Clean the skin with a disinfecting solution.
- A2.** Lay the drape with opening in position.

B. INFILTRATION OF LOCAL ANAESTHESIA

The local anaesthetic is NOT included in the SEFFIHAIR® kit. Using a local anaesthetic containing

adrenalin is recommended.

B1. Choose the spot for insertion of the cannula: it should be 1 cm out of the perimeter of the harvesting area (about 20 cm in diameter).

B2. Use the 27G needle to inject the anaesthetic solution.

B3. Make an opening in the skin with the 18G needle.

B4. Connect the multi-hole harvesting cannula with the 10 ml luer lock syringe containing local anaesthetic.

B5. Insert the syringe with the cannula into the guide.

B6. Rotate the syringe so that the barrel flange gets caught in the lock at the back of the guide.

B7. Make sure the cannula is parallel to the guide rod; if not so, set it in position with a delicate manual move.

B8. Insert the tip of the cannula into the skin opening perpendicularly until the guide stops.

B9. Rotate the cannula by 90 degrees so as to set the external guide parallel to the skin plane.

B10. Infiltrate the local anaesthesia with a back and forth movement of the cannula while maintaining the guide parallel and adhering to the skin. The anaesthetic solution should be injected in the whole harvesting area. If more anaesthetic is necessary, repeat the same procedure from step **B5** after recharging your syringe with more anaesthetic solution.

B11. WAIT 15 MINUTES.

C. HARVESTING OF TISSUE

C1. Connect the multi-hole harvesting cannula with the 10 ml luer lock syringe.

C2. Insert the syringe plunger lock at the back of the guide as shown in the pictures.

C3. Insert the syringe with the cannula into the guide.

C4. Rotate the syringe so that the barrel flange gets caught in the lock at the back of the guide.

C5. Make sure the cannula is parallel to the guide rod; if not so, set it in position with a delicate manual move.

C6. Insert the tip of the cannula into the skin opening perpendicularly until the guide stops.

C7. Rotate the cannula by 90 degrees so as to set the external guide parallel to the skin plane and the cannula in the subcutaneous plane.

C8. While holding the guide with your dominant hand, pull the syringe plunger back so as to create the vacuum in the syringe barrel until the plunger is caught in the back latch.

C9. While maintaining the guide parallel and adhering to the skin, move the cannula back and forth in the whole harvesting area. The tissue harvested will be sucked into the barrel and can be checked in quantity through the window in the upper body of the guide. Stop as soon as you have collected 5 ml of tissue in the syringe.

C10. Pull out the cannula from the skin.

C11. Rotate the syringe until the plunger flange is released from the lock at the back of the guide.

C12. Pull out the syringe from the guide and disconnect the cannula.

C13. Fill the syringe with sterile physiological solution (NOT included in the kit) you previously poured into the cylindrical container.

C14. Use the luer lock plug to close the syringe.

C15. Put the syringe (blue plug down) in the multiple syringe decanting stand.

Repeat the procedure from step **C1** to step **C15** until you have from one to three 10 ml syringes of tissue.

PREPARATION OF TISSUE

After some minutes the syringes in the decanting stand will show a separation by gravity of the tissue (top) from the physiological solution (bottom).

D. WASHING THE TISSUE

D1. Pull out the syringe from the multiple stand, hold it in a vertical position (DO NOT TILT NOR TURN UPSIDE DOWN), remove the plug, discard the washing physiological solution while retaining all the tissue inside the syringe.

D2. Put the plug back in position and the syringe back in the stand.

D3. Repeat the steps **D1** and **D2** for all the syringes in the stand.

E. FLUIDIFICATION

E1. Pull out the 10 ml syringe containing the tissue from the stand and remove the blue plug.

E2. Use the connector to connect the syringe to an empty 10 ml syringe.

E3. After pushing the tissue about 20 times from one syringe to the other, a good fluidification of the tissue is obtained.

E4. Disconnect the syringes and use the plug to close the syringe with the tissue.

E5. Put the syringe in the stand.

E6. Repeat the procedure from step **E1** to step **E5** for all the syringes.

WARNING:

The fluidification is to be considered a mechanical manipulation to make the tissue more fluid (and hence more easily implanted with a thinner needle), however it causes a varying loss of vitality in the adipose derived and mesenchymal cellular component.

F. CENTRIFUGATION

Centrifugation is performed when the therapy requires liquid tissue with no adipose stromal part. (Studies have demonstrated a small SVF reduction in the obtained liquid as compared to the whole tissue).

F1. Pull out the syringe from the multiple stand, hold it in a vertical position.

F2. Transfer the tissue from the syringe into the test tube.

F3. Close the test tube cap with the blue luer lock plug.

F4. Insert the test tubes into the centrifuge slots.

F5. Repeat the procedure from step **F1** to step **F4** for all the syringes in the stand while paying attention to setting the 3 test tubes in alternating positions inside the 6-slot centrifuge.

F6. Centrifuge the tissue for 10 minutes at 3000 RPM.

F7. Remove the test tubes from the centrifuge holding them in a vertical position.

G. TRANSFERRING TISSUE TO GRAFTING SYRINGES

G1. Remove the blue cap.

G2. With the long needle suck the liquid tissue from the bottom of the test tube, leaving the adipose stromal part and the oily component in the test tube.

G3. Use the blue luer lock plug to close the 3 ml syringe.

G4. Repeat the procedure from step **G1** to step **G3** until all the tissue contained in the test tubes is transferred to 3 ml syringes.

GRAFTING

H1. Remove the plug from the 3 ml syringe.

H2. Disinfect thoroughly the scalp.

H3. Insert the short mesotherapy needle (part of the kit).

H4. Inject the tissue in the scalp according to the treatment indications.

DISPOSAL

In case of use on patients with contagious bacterial and/or viral infections, some components could result highly contaminated.

Once used, do not dispose of the device with domestic waste. The device must be treated as a **Special Waste**.

Please deliver your used device to an authorised collection center.


Contact the Local Authorities to learn about the disposal regulations of such waste.

In order to respect the environment and not to cause any damage to health, we ask

you to comply with these recommendations.

The presence of this symbol highlights the obligation of differentiated waste disposal.

SEFFIHAIR®

The new age of
REGENERATIVE MEDICINE FOR HAIR LOSS THERAPY

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The product SEFFILLER® complies
with the requirements of the directive
93/42/CEE on medical devices.

SEFFILINE

**Via delle Lame, 98
40122 Bologna - Italy
+39 051 6494945**