

MRI Compatibility

ATTESTATION

I undersigned, Franck PELLETIER, certify that the Adjustable Gastric Bandings BIORING[®], Slimband[™] and ADHESIX[®] BIORING[®] manufactured by Cousin Biotech are compatible and safe with Magnetic Resonance Imaging (MRI).

Indeed the implanted raw materials used in the Adjustable Gastric Bandings BIORING[®], Slimband[™] and ADHESIX[®] BIORING[®] are non ferromagnetic materials.

Brand Name of adjustable gastric bandings	References	Implanted raw materials	MRI Compatibility
BIORING [®]	RINGS20000	Dimethyl siloxane + Polyethylene terephthalate + Barium sulphate + Polyoxymethylene	YES
	BCBRI23POU		YES
	BCBRIXLPOU		YES
Slimband [™]	BSBRI23POU		YES
	BSBRIXLPOU		YES
ADHESIX [®] BIORING [®]	BCBRI20ADU		Dimethyl siloxane + Polyethylene terephthalate + Barium sulphate + Polyoxymethylene + Polypropylene + Polyethylene Glycol + Polyvinylpyrrolidone
	BCBRI23ADU	YES	
	BCBRIXLADU	YES	
BIORING [®]	RINGPORTP1	Polysulfone + Titanium + Dimethyl siloxane + Barium sulphate	YES



Mr. Franck PELLETIER
Regulatory Affairs Director
January 08th, 2015

This document is confidential and the exclusive property of the firm COUSIN BIOTECH and cannot be reproduced, communicated or executed without its written permission

Headquarters : 8, rue de l'Abbé Bonpain - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 19 - **Fax** +33 (0) 3 20 14 40 13
 Factory : Allée des roses - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 20 - **Fax** +33 (0) 3 20 14 40 45
www.cousin-biotech.com - **e-mail** f.pelletier@cousin-biotech.com